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
Solid autosampler FPG 48

Automatic TIC Solids Module

Fully automated TOC determination with multi EA® 4000 thanks to TIC Solids Module “automatic”

TOC (Total Organic Carbon) and TIC (Total Inorganic Carbon) analyses are faster and easier than before! Easily determine the TOC and TIC parameters from one single sample in only one analysis step.

In compliance with standard methods: DIN EN 13137, prEN 15936, ISO 10694.

How does it work?

The TIC Solids Module “automatic” enables the automatic determination of Total Inorganic Carbon (TIC) in solid samples in combination with the sampler FPG 48. This way it is also possible to automate determination of Total Organic Carbon (TOC) either by difference or direct method. Additional sample pre-treatment like acidification is now redundant!

What can be determined?

It is possible to run 3 different method. All in one cycle – weigh in sample - Done!

TIC:

The sample is weighed into a single ceramic boat and placed on the sampler rack. The boat is automatically moved forward to the acidification position (outside of the furnace) where the acid is added drop-wise to the sample. The formed CO₂ is transferred to the NDIR detector. After the reaction is completed the ceramic boat is brought back to the sampler rack or alternatively to a waste collecting basin.

TOC by difference method:

It is the most preferred and recommended method for common applications. The same sample is weighed in two portions into two ceramic boats. The first boat is acidified like described above – TIC determination is performed. The second boat is forwarded completely into the furnace where combustion/oxidation takes places. This way the TC value is obtained. After the combustion is completed, the ceramic boat is brought back to the sampler rack or alternatively to a waste collecting basin.

In conclusion the TOC value is calculated automatically by the software as difference of TC-TIC.

TOC by direct method:

The sample is weighed into a single ceramic boat and placed on the sampler rack. The boat is automatically moved forward to the acidification position (outside of the furnace) where the acid is added drop-wise to the sample. The formed CO₂ is transferred to the NDIR detector.

After the reaction is completed the boat is moved forward into the furnace where combustion/oxidation takes places. This way the TOC value is directly obtained.

Solid autosampler FPG 48